

Package ‘metaR’

December 5, 2017

Type Package

Title Read Micromanager Metadata

Version 0.1.0

Author Julien Godet

Maintainer Julien Godet <julien.godet@unistra.fr>

Description Micromanager generates metadata files storing microscopes acquisition parameters. metaR gathered a collection of functions to easily read and extract these data.

License GPL-3

Encoding UTF-8

LazyData true

RoxygenNote 6.0.1

R topics documented:

aduToPhotons	2
getAcqTime	2
getDcamExposure	3
getElapsedTime	3
getEmGain	4
getExposure	4
getFrameNumber	5
getImageHeight	5
getImageWidth	6
getMetaGeneric	6
getParamNames	7
getXpos	7
getYpos	8
getZpos	8
Index	9

aduToPhotons	<i>Return photons from ADU</i>
--------------	--------------------------------

Description

Compute signal in photons from camera ADU. Set for Hamamatsu cam ImagEM.

Usage

```
aduToPhotons(EMGain)
```

Arguments

EMGain	EM gain from the camera
--------	-------------------------

Author(s)

JuG

getAcqTime	<i>Get Acquisition time - POSIXct object</i>
------------	--

Description

Get Acquisition time

Usage

```
getAcqTime(path)
```

Arguments

path	Path to the .txt or .gz metadata file
------	---------------------------------------

Value

Acquisition time (POSIXct object)

Author(s)

JuG

getDcamExposure	<i>Get DCAM Exposure</i>
-----------------	--------------------------

Description

Get DCAM Exposure - expressed in msec

Usage

```
getDcamExposure(path)
```

Arguments

path	Path to the .txt or .gz metadata file
------	---------------------------------------

Value

DCAM exposure in ms

Author(s)

JuG

getElapsedTime	<i>Function to extract elapsed time (in msec)</i>
----------------	---

Description

Get elapsed time in ms between consecutive frames

Usage

```
getElapsedTime(path)
```

Arguments

path	Path to the .txt or .gz metadata file
------	---------------------------------------

Value

Elapsed time (vector) in msec

Author(s)

JuG

`getEmGain`*Get EM Gain*

Description

Get EM Gain

Usage`getEmGain(path)`**Arguments**

<code>path</code>	Path to the .txt or .gz metadata file
-------------------	---------------------------------------

Value

EMGain value

Author(s)

JuG

`getExposure`*Get Exposure time in ms*

Description

Get Exposure time in ms

Usage`getExposure(path)`**Arguments**

<code>path</code>	Path to the .txt or .gz metadata file
-------------------	---------------------------------------

Value

Exposure time in ms

Author(s)

JuG

getFrameNumber	<i>Get the total number of frames</i>
----------------	---------------------------------------

Description

Get frame number

Usage

getFrameNumber(path)

Arguments

path	Path to the .txt or .gz metadata file
------	---------------------------------------

Value

EMGain value

Author(s)

JuG

getImageHeight	<i>Get Image Height</i>
----------------	-------------------------

Description

Get image height - expressed in pixels

Usage

getImageHeight(path)

Arguments

path	Path to the .txt or .gz metadata file
------	---------------------------------------

Value

Image height in px

Author(s)

JuG

getImageWidth	<i>Get Image Width</i>
---------------	------------------------

Description

Get image width - expressed in pixels

Usage

```
getImageWidth(path)
```

Arguments

path	Path to the .txt or .gz metadata file
------	---------------------------------------

Value

Image width in px

Author(s)

JuG

getMetaGeneric	<i>Generic function to get a given parameter from the metadata file</i>
----------------	---

Description

Generic function to get a given parameter from the metadata file need to define time zone using Sys.setenv(TZ='GMT') - otherwise warning message

Usage

```
getMetaGeneric(path, parameter = "FrameIndex", type = "ai")
```

Arguments

path	Path to the .txt or .gz metadata file
parameter	parameter to extract (string)
type	i:integer, n:numeric, f:factor, t:time, ai:'asis'

Value

Values of the parameter for each image

Author(s)

JuG

Examples

```
pNames <- getParamNames(path)
pNames[59]
getMetaGeneric(path, parameter = pNames[59], type = "ai")

getMetaGeneric(path, parameter = "Time", type = "t")
```

getParamNames*Get Metadata parameters names*

Description

Get metadata parameters names

Usage

```
getParamNames(path)
```

Arguments

path Path to the .txt or .gz metadata file

Value

Metadata parameters names

Author(s)

JuG

getXpos*Function to extract X position in micrometer*

Description

Get X position in micrometer

Usage

```
getXpos(path)
```

Arguments

path Path to the .txt or .gz metadata file

Value

X position (vector) of images

Author(s)

JuG

`getYpos`*Function to extract Y position in micrometer*

Description

Get Y position in micrometer

Usage

```
getYpos(path)
```

Arguments

path Path to the .txt or .gz metadata file

Value

Y position (vector) of images

Author(s)

JuG

`getZpos`*Get Z position in micrometer*

Description

Get Z position in micrometer

Usage

```
getZpos(path)
```

Arguments

path Path to the .txt or .gz metadata file

Value

Z position (vector) of images

Author(s)

JuG

Index

aduToPhotons, [2](#)

getAcqTime, [2](#)

getDcamExposure, [3](#)

getElapsedTime, [3](#)

getEmGain, [4](#)

getExposure, [4](#)

getFrameNumber, [5](#)

getImageHeight, [5](#)

getImageWidth, [6](#)

getMetaGeneric, [6](#)

getParamNames, [7](#)

getXpos, [7](#)

getYpos, [8](#)

getZpos, [8](#)